
Green RUSCH® Lite™

Single-Use Metal Laryngoscope Blade for Fiber Optic System

en Single-Use Metal Laryngoscope Blade for Fiber Optic System READY FOR USE.

INSTRUCTIONS FOR USE: 1. Attach to Fiber Optic system handle. 2. Click into place. 3. Lift blade and verify illumination. 4. After use, dispose according to standard biological hazard procedures.

INTENDED USE: The Laryngoscope is a device intended to allow direct visualization of the laryngeal structure and is used during tracheal intubation.

WARNING: Single use: Do not reuse or reprocess. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

SIDE EFFECTS:

- Laryngeal/airway trauma
- Bleeding
- Sore throat
- Hoarseness
- Dental trauma
- Dysphagia

Do not sterilize or autoclave. This product is not suitable for an MRI environment. Do not use if individual unit has been opened or damaged. Only personnel adequately trained to perform tracheal intubation should use laryngoscopes.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

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