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Item: AEROTube® – Disposable PVC Laryngeal Mask

1. Product Description

The AEROTube® – Laryngeal Mask is inserted into the trachea through glottis to establish and maintain an artificial airway during anaesthesia and emergency treatment.

2. Intended Users

Nursing staff, medical professionals.

3. Indications

The product is intended for use in the emergency department, ICU, or other departments for anesthesia, first aid or resuscitation etc. It is also possible to use it in situations where manipulation of the head or neck to facilitate endotracheal intubation is difficult or unsuccessful.




4. Contraindication

The AEROTube® – Laryngeal Mask may only be used on patients who, according to medical opinion, are suitable for application of AEROTube® – Laryngeal Mask.

The use of AEROTube® – Laryngeal Masks is strictly prohibited in the following cases / situations:

- 4.1 Patients at increased risk of aspiration or regurgitation and non-fasted patients.
- 4.2 Excessive abdominal pressure, vomit and reflux, habitual vomiting.
- 4.3 Throat infection or other pathological changes.
- 4.4 Bleeding in the respiratory tract.
- 4.5 Tonsil abnormal enlargement, potential airway obstruction, throat sore, abscess, and hematoma is strictly prohibited.
- 4.6 Patients, who are more than 14 weeks pregnant.
- 4.7 Patients with a firm, reduced lung compliance.
- 4.8 Patients who are extremely or pathologically overweight.
- 4.9 Patients where the mouth can not be opened sufficiently to allow insertion of the AEROTube® – Laryngeal Mask.
- 4.10 Not deeply unconscious patients, who resist the insertion of the AEROTube® – Laryngeal Mask.

5. Application Description

- 5.1 Choose the right tube size according to the patient's age, gender, height, weight etc.
- 5.2 Deflate the cuff completely so that it forms smooth "spoon-shape". Lubricate the posterior surface of the mask with water-soluble lubricant.
- 5.3 Hold the AERObute® – Laryngeal Mask like a pen, with the index finger placed at the junction of the cuff and the tube.
- 5.4 With the head extended and the neck flexed, carefully flatten the AERObute® – Laryngeal Mask tip against the hard palate.
- 5.5 Use the index finger to push cranially, maintaining pressure on the tube with the finger.
- 5.6 Advance the mask until definite resistance is felt at the base of the hypopharynx.
- 5.7 Gently maintain cranial pressure with the non-dominant hand while removing the index finger.
- 5.8 Without holding the tube, inflate the cuff with just enough air to obtain a seal (to a pressure of approximately 60 cm H₂O). See the instructions for appropriate volumes. Never overinflate the cuff.

Size	Patient Selection Information	Max. Inflation Volume (ml)
1	Neonates up to 5 kg	4
1.5	Infants 5-10 kg	7
2	Infants / Children 10-20 kg	10
2.5	Children 20-30 kg	14
3	Children 30-50 kg	20
4	Adults 50-70 kg	30
5	Adults 70-100 kg	40

6. Warning Notices

- 6.1 Be sure all connections are secure.
- 6.2 Deflate the cuff completely during the intubation and extubation.
- 6.3 Keep the tube airway unobstructed and timely aspirate secretions.
- 6.4 The patients with tracheal intubation should receive oral cleaning and care regularly.
- 6.5 Strengthen the management of tube airway's temperature and humidity.
- 6.6 During the indwelling time, release and re-inflate the cuff regularly to alleviate the pressure of trachea for a long time.
- 6.7 Tracheal intubation is intended for use in a short time. If further treatment is needed, other forms of intubation are recommended.

7. LOT Number

The LOT number is printed on the packaging.

8. Storage Conditions

Do not store in direct sunlight, at extreme temperatures or humidity. Store in dry conditions, in a non-corrosive atmosphere and in a well-ventilated room.

9. Handling Conditions

In case of damaged packaging or inappropriate storage the product **SHOULD NOT** be used and has to be destroyed.

10. Sterile Condition

The product is **STERILE**.

11. Sterilisation Process

Single use only. **DO NOT** re-sterilize.

12. Single Use Notice

The product is a disposable product.

13. Notice on Reuse

DO NOT reuse.

14. Reprocessing Notice

The product **SHOULD NOT** be reprocessed.


15. Notice on the Application by Laymen

The product is intended for use by qualified medical personnel.

16. Notice on Reporting Requirements

Please inform the manufacturer, the distributor and the importer in all cases of serious incidents in connection with the product.

17. Symbols

	Manufacturer		Mag. Date
	REF No.		EXP. Date
	LOT No.		Instruction
	EO sterilization		Latex free
	Don't reuse		Warning
	Don't use package damaged		Mask size
	Keep dry		Keep away from sunlight
	Do not re-sterile		Medical Device
	EC REP		European importer
	CE mark		Temperature limit
	Humidity, limit		Contains phthalates

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Hangzhou Shanyou Medical Equipment Co., Ltd.
Prolinx GmbH, ADD.: Brehmstr. 56,