

+33 (0)4.68.68.39.68 Fax: +33 (0)4.68.68.39.69

info@aivia.eu

WALL MOUNTED AIVIA PRE-INSTALLATION FILE






- Never install Aivia under direct sunlight exposure. You risk exposing the defibrillator to excessive temperatures.
- The Aivia is « Safety Extra-low voltage », for human security (The Aivia 100 is powered by 4 D-Cell Alkaline batteries).
- Setup and maintenance must comply by existing regulations.
- Pyrescom SAS releases it's entire responsability in an event of non-compliance with the recommandations.

A-DPM

V2.05 1/2 2019-02-11



Aivia installation must be done by qualified personnel only, and member of the Bornavie charter.

Ensure reliability of the support. Fixing must be adapted to the type of support the Aivia is being fixed on.

- Aivia implementation zone must be free as signified above.
- If the Aivia is installed in the street or any unrestricted passage way, envisage a ground surface (abutment) or a column.
- Drilling template is part of the Aivia packing. It is also available in the document named A-DPM-GABARIT.
- Aivia fixing must be able to carry a minimal load of 20 kg (44 lb).

POWER SUPPLY UNIT:

It must be outside the Aivia, never install it into the Aivia.

The power supply unit is made up from:

- One earth bar.
- 10A/30mA Circuit breaker.
- 24VDC/3A CLASS II Power supply*.
- The power supply must be limited in output power.
- Its mains wiring must contain Live, Neutral, and Earth/Ground.
- Full load voltage drop on the secondary must not exceed 2%, in that case 23,5V at the Aivia connector
- Power supply unit parts must be protected against water, dust, and ventilated**.
- Wiring of the power supply unit must comply by existing regulations in the country of installation.

SELV WIRING:

- U1000R2V Type. 2 Isolated and gained conductors, 1,5mm² gauge, 10 meters (393 in) maximum length.
- For +24V use a red or brown cable.
- For 0V use a black or grey cable.
- For earth connecting, use a green/yellow cable. Only STN connected Aivia must be connected to earth.

Recommanded power supply: PYRESCOM C AL/HDR100-24CA

Natural or forced convection en depending on the environnment and thermal caracteristics at maximum loa

Wiring: 4 pairs phone cable.

PSTN Analog*: Direct, or on private branch exchange.

(Aivia 220 230 320 330 °

GSM*: Ensure of the coverage area by the designated mobile network operator on the installation site.

Subscribtion being at the client's expense

INTERNET ACCES

- Ethernet network*: having an internet connectivity. HTTP protocol port: 80. Wiring: 1 Ethernet cable, CAT5 FTP or better.
- **GPRS*:** Ensure of the coverage area by the designated mobile network operator on the installation site.

Subscribtion being at the client's expense.