i-view™  
Video Laryngoscope






**Instructions For Use**

# **i-view™**

## **Video Laryngoscope**

<b>en</b>	<b>Instructions for use</b>
<b>fr</b>	<b>Mode d'emploi</b>
<b>de</b>	<b>Gebrauchsanweisung</b>
<b>es</b>	<b>Instrucciones de uso</b>
<b>pt</b>	<b>Instruções de utilização</b>
<b>it</b>	<b>Istruzioni per l'uso</b>
<b>nl</b>	<b>Gebruiksaanwijzing</b>
<b>no</b>	<b>Instruksjoner for bruk</b>
<b>fi</b>	<b>Käyttäjän opas</b>
<b>sv</b>	<b>Bruksanvisning</b>
<b>da</b>	<b>Brugsanvisning</b>
<b>el</b>	<b>Οδηγίες χρήσης</b>
<b>lt</b>	<b>Naudojimo vadovas</b>
<b>pl</b>	<b>Instrukcja użytkowania</b>
<b>ru</b>	<b>Инструкция по применению</b>
<b>cs</b>	<b>Návod k použití</b>
<b>hu</b>	<b>Használati utasítás</b>
<b>sl</b>	<b>Navodila za uporabo</b>
<b>lv</b>	<b>Lietošanas instrukcija</b>
<b>et</b>	<b>Kasutamjuhised</b>
<b>bg</b>	<b>Инструкции за употреба</b>
<b>sr</b>	<b>Uputstvo za upotrebu</b>
<b>ro</b>	<b>Instrucțiuni de utilizare</b>
<b>sk</b>	<b>Návod na použitie</b>
<b>hr</b>	<b>Upute za upotrebu</b>
<b>tr</b>	<b>Kullanma kılavuzu</b>
<b>ja</b>	<b>取扱説明書</b>
<b>zh-s</b>	<b>使用说明</b>
<b>zh-t</b>	<b>使用說明</b>
<b>ar</b>	<b>تعليمات الاستخدام</b>

The medical device is not supplied individually. Only one copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users.

**NOTE:** Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product.

**WARNING:** A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION:** A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

### Description

- The Intersurgical i-view™ is a single use, hand-held fully disposable video laryngoscope, suitable for use with adult patients to facilitate oral tracheal intubation, including during anaesthesia and emergency medicine.

### Intended use

To facilitate laryngoscopy.

### Warnings

- For use by qualified medical personnel only. Operators must familiarise themselves with these instructions for use, and the device.
- No modification of this device is permitted.
- The device contains no user serviceable parts inside.
- The batteries within the device should not be removed prior to its disposal.
- No attempt should be made to reprocess, sterilise or re-use the device.
- Do not use force to insert the device.
- Particular care should be taken with patients who have fragile and vulnerable dental work, in accordance with recognised airway management practices and techniques.
- In circumstances where fogging is regarded as a possibility, such as when the patient is exposed to high flow humidified oxygen before or during use of the i-view, the device should be turned on > 30 seconds prior to use.
- To reduce the potential for trauma, the i-view should always be inserted into the patient's mouth under direct vision.
- To reduce the potential for trauma, the endotracheal tube should always be initially introduced into the patient's mouth under direct vision. The endotracheal tube may then be advanced through the vocal cords by direct or indirect vision.
- Ensure that a backup plan is in place in case of difficulty or emergency while performing the procedure.
- Do not use the device in strong magnetic fields such as MRI environments.
- Users should use both direct and indirect (video) visualisation to progressively expose the key anatomical landmarks in accordance with best practice for laryngoscopy.
- When turning off the device, the on/off button should be depressed for more than half a second.
- During prolonged use, the clear lens in front of the light source can reach temperatures up to a maximum of 50°C and should not be held in direct contact with patient tissue for longer than 1 minute.

### Indications For Use

i-view is designed to provide direct and indirect visualisation of the larynx. It helps to facilitate oral tracheal intubation in adults by personnel who are suitably trained and experienced in the use of airway management techniques and devices, including intubation. i-view is approved for use in road ambulances under IEC 60601-1-12.

### Contraindications

- Surgical airway.
- Fluid or solid material in the airway that cannot be cleared with suction.
- Trismus, limited mouth opening, pharyngo-perilaryngeal abscess, trauma or mass.
- Operator inexperience with the device.

### Size Selection

- i-view is available in one size for adult patients. In smaller adult patients, it may not be necessary to insert the blade as far as may be necessary for an average sized adult patient in order to gain an optimal view of the relevant anatomical structures. For patients where the blade of the i-view is not large enough to enable these structures to be viewed, an alternative device should be used.

### Pre-use checks

- Inspect the packaging and ensure it is not damaged prior to opening.

- Inspect the device carefully and discard if the device looks abnormal or deformed.
- Turn on the device by depressing the on/off switch (**fig 1**) and check there is a visible light at the distal end of the blade and that the Intersurgical logo appears on the screen for a few seconds. Turn off the device by depressing the on/off button. The device is now fully operational. Discard the device if it cannot be successfully switched on and off.
- In order to reduce the potential for battery depletion, i-view will automatically switch itself off after 20 minutes. If required, the device can be switched on again by pressing the on/off button.

#### Using the i-view

- Turn on the device by pressing the on/off button (**fig 1**).
- The distal tip of the blade of the device should be inserted into the patient's mouth under direct vision (**fig 2**) and the blade gently advanced down the tongue until the epiglottis has been identified (**fig 3**). The tip of the blade should then be advanced into the vallecula (**fig 4**).
- With minimal force, indirectly lift the epiglottis with the i-view until the glottis comes into view. This can be achieved using a direct or indirect view (**fig 5**). If deemed necessary, bimanual laryngoscopy and/or head elevation can be used to improve the view.
- Once the glottis has been identified, an appropriately sized tracheal tube can be passed through the vocal cords to the required depth (**fig 6**) and the i-view gently removed (**fig 7**).
- Confirmation of correct placement of the tracheal tube should be assessed by the use of wave form capnography and other methods in accordance with applicable guidelines.
- i-view should then be switched off by depressing the on/off switch (**fig 8**).

#### Storage Conditions

Store in a dry place out of direct sunlight. Products should be stored between -18 °C - +40 °C.

#### Disposal

- Should be disposed of in accordance with normal procedures for clinical waste.
- Please contact your local authority or Intersurgical for more information.

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