Arrow® T-POD™ Pelvic Stabilization Device

Rx only.

Indications for Use:

Provides circumferential compression to the pelvis in patients with suspected pelvic fracture for pelvic stabilization which may reduce blood loss and pain.

Disposal:



Where applicable dispose of according to local regulations for biohazardous waste.

When Used In X-Ray:

Detachable X-Ray Detectable Label is included to identify that T-POD is in place on patient.

⚠️ General Warnings

1. Warning: Single use: Do not reuse, reprocess or sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Photos are for illustration purposes.

Instruction For Use:

1. Slide the Belt under the patient and into position under the pelvis. It is recommended to remove patient clothing prior to applying T-POD.
2. Trim the Belt leaving a 6-8" (15-20 cm) gap over the center of the abdomen. Note: If unable to trim, the Belt can be folded under.
3. Apply the Velcro®-backed Pulley System to the Belt on each side of the gap.
4. Draw the Pull Tab, creating simultaneous circumferential compression.
5. Wrap the cord around the Pulley System hooks and secure the Velcro®-backed Pull Tab to the Pulley System fabric.
6. Record the date and time of application.

Symbol Glossary

								
Medical device	Caution	Consult instructions for use	Do not reuse	Non-Sterile	Keep away from sunlight	Keep dry	Do not use if package is damaged	Not made with natural rubber latex
								
MR Safe	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture			

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

Teleflex, the Teleflex logo, Arrow, the Arrow logo and T-POD are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Other names may be trademarks of their respective owners. © 2019 Teleflex Incorporated. All rights reserved.

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