- de** Gebrauchsanleitung
FEROtube® – Endotrachealtuben
- fr** Mode d'emploi
FEROtube® – Tubes endotrachéaux
- it** Istruzioni per l'uso
FEROtube® – Tubi endotracheali
- pt** Instruções de utilização
Tubos endotraqueais – FEROtube®
- ro** Instrucțiuni de utilizare
FEROtube® – Tuburi endotraheale
- sk** Návod na použitie
FEROtube® – Endotracheálne hadičky
- el** Οδηγίες χρήσης
Ενδοτραχειακοί σωλήνες – FEROtube®
- sl** Navodila za uporabo
FEROtube® – Endotrahealni tubusi
- sv** Bruksanvisning
FEROtube® – Endotrakealtuber
- lt** Naudojimo nurodymai
FEROtube® – Endotrachėjiniai vamzdeliai
- mt** Istruzzjonijiet għall-użu
Pajpijiet Endotrakeali – FEROtube®
- hu** Használati utasítás
FEROtube® – Endotracheális tubusok

- en** Instruction for use
FEROtube® – Endotracheal Tubes
- nl** Gebruiksaanwijzing
FEROtube® – Endotracheale buizen
- es** Instrucciones de uso
Tubos endotraqueales – FEROtube®
- pl** Instrukcja użytkowania
FEROtube® – Rurki dotchawicze
- cs** Návod k použití
FEROtube® – Endotracheální hadičky
- bg** Инструкции за употреба
FEROtube® – ендотрахеални тръби
- hr** Upute za upotrebu
FEROtube® – Endotrahealni tubusi
- da** Brugsanvisning
FEROtube® – Endotrakealtuber
- fi** Käyttöohjeet
FEROtube® – Endotrakeaaliputket
- lv** Lietošanas pamācība
FEROtube® – Endotracheālās caurules
- ga** Treoir maidir le húsáid
FEROtube® – Tiúba Iontraicéacha
- et** Kasutusjuhend
Endoträhhaaltorud – FEROtube®

FEROtube® – ET01

FEROtube® – ET01-M

FEROtube® – ET01-M-SET

de Inhalt

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en English

REF: HET01-... | HET01-M-... | HET01-M-SET-...

Item: AEROTube® – Endotracheal tube, uncuffed | AEROTube® – Endotracheal tube, cuffed |
AEROTube® – Endotracheal tube-set, cuffed

1. Models

This instruction for use applies to different types of AEROTube® – endotracheal tubes.

- Standard AEROTube® – endotracheal tubes:
 - Without cuff
 - With cuff
 - With cuff, stylet and Luer syringe
- Reinforced AEROTube® – endotracheal tubes:
 - Without cuff
 - With cuff
 - With cuff, stylet and Luer syringe

Not all versions are available in all models and sizes.

2. Intended Use

- Remove the sterile tube from the protective packaging.
- If necessary, the tube may be shortened prior to intubation by removing the tube slid halfway onto the connector, cutting it to the correct length and sliding the connector back firmly into the tube. Please follow the safety precautions.
- Ensure that the tube is correctly connected with the ventilation system.
- Block-test the cuff, control cuff and valve for airtightness before use. To do this, fill a syringe with a conical Luer taper with air, insert it firmly into the valve until an airtight connection is established and fill the cuff until the cuff skin is taut.
- Then empty the cuff completely.
- Following intubation, fill the cuff via the Luer syringe or using an appropriate cuff pressure gauge until an effective seal is obtained. Please follow the safety precautions.
- Remove the syringe from the valve once the cuff is filled, as otherwise the valve will remain open and air may leak from the filled cuff.
- Ensure that the system is absolutely airtight. Monitor cuff pressure carefully. If there are any changes in pressure, the cause must be determined and remedied immediately.
- Prior to extubation, all air must be removed from the cuff via the Luer syringe.
- AEROTube® – endotracheal tubes must be disposed of in accordance with national guidelines for hospital hygiene and waste disposal.

3. Indications

- AEROTube® – endotracheal tubes are used to assist respiration and supply air to the respiratory system during anesthesia, intensive care and emergency medical care.
- They are introduced via endotracheal intubation and are thus used for short- and long-term nasal/oral intubation.
- They can also be used to administer drugs nasally and/or orally.

4. Contraindications

- None known.

5. Description

- AERObute® – endotracheal tubes are manufactured from medical-grade polyvinyl chloride (PVC).
- The tube is slid only halfway onto the connector to enable the length of the tube to be adjusted to fit the patient.
- Different models of tubes are available with and without a cuff and with and without a metal spiral embedded.
- Using a cuff is intended to prevent any blood, secretions or foreign bodies entering the airways and thus being aspirated.
- The filling system for the cuff consists of an inlet tube, a control cuff (blue) and a self-closing valve with a conical Luer taper.
- An inbuilt lateral opening under the cuff (Murphy eye) ensures that sufficient resuscitation can be performed even if the end of the tube is blocked.
- AERObute® – endotracheal tube sets include a guide pin (stylet) in the tube and a Luer syringe (10 ml) as part of the set.
- The integral guide pin (stylet) is intended to facilitate the intubation of the tracheal tube and the Luer syringe to fill the cuff.
- Reinforced endotracheal tubes may be used to reduce the potential for kinking in case of unusual positioning of the head and neck following the intubation.

6. Safety instructions for endotracheal tubes

- If a laser or surgical probe is used during treatment, there is a risk, particularly where oxygen or nitrous oxide/air mixtures are used, that the tube will ignite, causing burns and releasing toxic, irritant vapors and gases, including hydrochloric acid (HCL).
- If the tube is shortened, the connector must not be placed in direct proximity of the input to the cuff filling system in the tube shaft, as this can render the device no longer airtight. Reinforced endotracheal tubes cannot be cut to length.
- The intubation markings (black rings or dots) help with the placement of the tube but, owing to differences in the configuration and length of airways from patient to patient, are not under any circumstances a substitute for a clinical assessment by a doctor.
- If the position of the patient or the tube changes following intubation, ensure that the tube is still in the correct position.
- The correct size of tube for each patient must always be determined in accordance with a clinical assessment.
- Intubation and extubation should be performed in accordance with current medical standards.
- If a lubricant is used, it is essential to ensure that the lumen is not blocked, which would prevent respiration and damage the cuff.
- Lubricant must not be used to facilitate any reinsertion of the connector, as this would increase the risk of the connector accidentally loosening/disconnecting.

7. Additional safety instructions for AERObute® – endotracheal tubes with cuffs

- Regardless of the type of cuff used, the user must ensure that the cuff is only ever filled with the minimum air volume necessary to achieve an effective seal at the desired respiratory pressure.
- Never fill the cuff without measuring the volume or with a predefined volume of air, as the pressure resistance in the control cuff alone is not a sufficient indicator of the pressure of the cuff.
- Never overfill the cuff! Cuff pressure should not as a general rule exceed 25 cm H₂O. Excess pressure can lead to damage to the trachea, blockage of the airways or damage to the cuff (which in turn will result in a loss of pressure).
- Diffusing nitrous oxide or oxygen can also change the filling volume.
- Using a cuff pressure gauge can assist in monitoring and regulating cuff pressure.
- Aerosol lubricants can cause damage to the cuff. The use of a water-soluble lubricant gel is recommended.
- Care should be taken to prevent damage to the cuff during intubation, e.g. from teeth, nasal conchae or any rough surfaces on the intubation instruments. AERObute® – endotracheal tubes with damaged cuffs must be removed from use.
- Unlock the cuff before making any change to the position of the tube, otherwise the locked cuff may injure the patient or be damaged.
- Auxiliary instruments such as syringes, shut-off valves, pressure gauges etc. should not be left in the cuff filling system for extended periods.

General warnings and precautionary measures

- AERObute® – endotracheal tubes must be used by trained and authorized specialist personnel only, who must have a sufficient level of training in their use.

- Intubation and extubation should be performed in accordance with current medical standards.
- The appropriate diameter of tube for each patient must be determined individually under the clinical supervision of a doctor.
- If intubation is performed using a guide pin (stylet), care must be taken to ensure that the tip of the stylet does not protrude beyond the end of the tube or through the Murphy eye.
- When using the AEROTube® – endotracheal tube sets, remove the guide pin (stylet) from the tube before intubation and moisten it with a suitable lubricant.
- If intubation is performed using a guide pin (stylet), moisten the stylet with a suitable lubricant before intubation.
- To protect the tube against damage from the patient's bite reflex, a bite block or an oropharyngeal tube should be put in place immediately after intubation.
- AEROTube® – endotracheal tubes contain phthalates – see product label.
- Phthalates can have potentially toxic effects on reproduction. According to current scientific knowledge, it cannot be excluded that long-term exposure or use may pose a risk to premature male babies. Medical products containing phthalates should be used only for short periods on pregnant women, breastfeeding mothers, babies and toddlers.

8. Storage instructions

- Optimum storage temperature: - 10° to + 40°C
- Keep away from sunlight and keep dry.
- Do not use if the product sterilization barrier or the packaging is compromised.
- Exposure to elevated temperature and UV light and excessive moisture should be avoided during storage.

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Co., Ltd., No. 18, Wenzhou Road, Economic
SUNGO Europe B.V., Fascinatio