## Quicktrach I

---

CE 0123

# Deutsch

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## VISUAL INSPECTION

- ▶ Check the product for damage and loose particles. A faulty product must be disposed of (see "Disposal").

## USE

- ▶ Extend the neck of the patient.
  - ▶ Locate the cricothyroid membrane (b) by the palpable depression between the thyroid cartilage (a) and cricoid cartilage (c).
  - ▶ Stabilise this site with index finger and thumb.
  - ▶ Puncture the cricothyroid membrane (b) perpendicular to the trachea.

Due to the conical needle tip, the opening into the trachea is dilated. An initial skin incision with a scalpel could reduce the force required.

### CAUTION

The posterior tracheal wall can be injured by inserting the Quicktrach I too deeply or by removing the stopper too early.

- ▶ Insert the Quicktrach I in caudal direction to the stopper.
  - ▶ Aspirate air with the syringe to determine the position of the cannula. A syringe pre-filled with e.g. normal saline (2 - 3 ml) can be helpful. If aspiration of air is possible, the needle tip is in the trachea.
- ▶ Remove the stopper from the cannula.

### CAUTION

The needle serves as a guide and must not be inserted further.

- ▶ While holding the needle in position, advance the cannula along the needle until the flange is flush with the patient's neck.
  - ▶ Remove the needle, while holding the cannula in position.
  - ▶ Ventilate the patient directly via the 15 mm connector or via the included extension tubing for enhanced flexibility.
- ▶ Secure the cannula with the necktape.

Quicktrach I provides a temporary airway and can be used to ventilate the patient's lungs until the airway can be secured by other techniques.

### Removal after use

- ▶ After use, remove the necktape and the cannula.

## SHELF LIFE

Use-by date: See product label.

## STORAGE AND TRANSPORT CONDITIONS

### CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

## DISPOSAL

### CAUTION

- The product may be contaminated with potentially infectious substances of human origin.
- The needle tip is sharp.

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

## PRODUCT SPECIFICATIONS

<b>REF</b>	30-04-002-1	30-04-902-1	30-04-004-1	30-04-904-1
<b>Patient</b>	Child		Adult	
<b>Needle</b>	2.0 mm I.D.		4.0 mm I.D.	
<b>Cannula connection</b>	15 mm O.D.			

## English

### INTENDED USE

Percutaneous cricothyrotomy set for securing the airway in case of upper airway obstruction.

Clinical benefit: ventilation of a patient

Patient target group: children and adults

Place of use: hospital and prehospital including military use

### INDICATION

Upper airway obstruction or ultima ratio if all other attempts to ventilate the patient have failed.

### CONTRAINDICATION

None known.

### SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them with the product.

- The product may only be used by a physician or authorised and trained medical personnel who have sufficient knowledge in handling the product.
- Only use the product for the patient target group.
- The needle tip is sharp.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the product must be inspected visually for damages (see "Visual inspection"). A faulty product must not be used.
- Maximum application time: 4 hours
- Constantly monitor ventilation through standard techniques or the rising and falling of the chest.
- It is not permitted to make any changes to the product.
- The product is not suitable for MRI.

- The product is intended for single use and must not be reused and/or reprocessed. The function of the product is impaired by reprocessing. Any reuse entails a potential infection hazard.
- The product is sterile (ethylene oxide).
- The product must not be used if the packaging is damaged or the expiration date has expired.

REF	30-04-002-1	30-04-902-1	30-04-004-1	30-04-904-1
Connection extension tubing to breathing system	15 mm O.D.	22 mm I.D.	15 mm O.D.	22 mm I.D.
Weight (incl. packaging)	140 g	140 g	145 g	145 g
Packaging unit	1 pc in blister pack	1 pc in plastic tube	1 pc in blister pack	1 pc in plastic tube

#### MATERIAL DATA

Cannula	ABS (acrylonitrile butadiene styrene), EVA (ethylene vinyl acetate copolymers)		
Stopper	ASA (acrylonitrile styrene acrylate)		
Needle	Stainless steel		
Syringe	PP (polypropylene)		
Necktape	Foam	Inside: PU (polyurethane), outside: PA (polyamide)	
	Hook and loop fastener	PES (polyester), PP (polypropylene)	
Extension tubing	30-04-002-1, 30-04-004-1: Silicone		
	30-04-902-1, 30-04-904-1: PP (polypropylene)		

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