
natural airway management

Instructions For Use

www.i-gel.com










Rx ONLY


natural airway management


i-gel® supraglottic airway	
en	Instructions for use
fr	Mode d'emploi
de	Gebrauchsanweisung
es	Instrucciones de uso
pt	Instruções de utilização
it	Istruzioni per l'uso
nl	Gebruiksaanwijzing
no	Brukerveiledning
fi	Käyttäjän opas
sv	Bruksanvisning
da	Brugsanvisning
el	Οδηγίες χρήσης
lt	Naudojimo vadovas
pl	Instrukcja użytkowania
ru	Инструкция по применению
cs	Návod k použití
hu	Használati utasítás
sl	Navodila za uporabo
lv	Lietošanas instrukcija
et	Kasutusjuhised
bg	Инструкции за употреба
sr	Uputstvo za upotrebu
ro	Instrucțiuni de utilizare
sk	Návod na použitie
hr	Upute za uporabu
tr	Kullanma kılavuzu
ja	取扱説明書
zh-s	使用说明
zh-t	使用說明
ar	تعليمات الاستخدام

Picture 1

Picture 2

The medical device is not supplied individually. Only one copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users.

NOTE: Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product.

WARNING: A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Product name & description

i-gel® supraglottic airway (adult sizes). Product codes: 8203000, 8204000 and 8205000.

Key components (picture 1)

- 1.1 Distal end of gastric channel
- 1.2 Epiglottic rest
- 1.3 Proximal end of gastric channel
- 1.4 Soft non-inflatable cuff
- 1.5 Buccal cavity stabiliser
- 1.6 Integral bite block
- 1.7 15 mm connector

Product positioning (picture 2)

Key:

- 2.1 Oesophagus
- 2.2 Trachea
- 2.3 Vocal cords/folds
- 2.4 Thyroid cartilage
- 2.5 Laryngeal Inlet
- 2.6 Sealing mechanism. Non-inflatable cuff
- 2.7 Ventilatory opening
- 2.8 Hyoid bone
- 2.9 Epiglottis
- 2.10 Tongue
- 2.11 Nominal depth of insertion mark
- 2.12 Incisor
- 2.13 Nominal length of the internal pathway (refer to technical data chart below)

Intended use

Securing and maintaining a patent airway. To provide a conduit for the delivery of respiratory gases to the patient.

Indications

The i-gel® is indicated for use in securing and maintaining a patent airway in routine and emergency anaesthetics of fasted patients, during spontaneous or intermittent positive pressure ventilation, during resuscitation of the unconscious patient, and as a conduit for intubation under fiberoptic guidance in a known difficult or unexpectedly difficult intubation, by personnel who are suitably trained and experienced in the use of airway management techniques and devices.

Contraindications

- Non-fasted patients for routine and emergency anaesthetic procedures.
- Trismus, limited mouth opening, pharyngo-perilaryngeal abscess, trauma or mass.
- Do not allow peak airway pressure of ventilation to exceed 40 cm H₂O.
- Do not use excessive force to insert the device or gastric tube.
- Inadequate levels of anaesthesia which may lead to coughing, bucking, excessive salivation, retching, laryngospasm or breath holding thus complicating the anaesthetic outcome.
- Do not leave the device in situ for more than 4 hours.
- Do not reuse or attempt to reprocess the i-gel®.
- Patients with any condition which may increase the risk of a full stomach e.g. hiatus hernia, sepsis, morbid obesity, pregnancy or a history of upper gastro-intestinal surgery etc.
- Use on a conscious/semi-conscious patient in an emergency setting.

Patient Category

The product range is suitable for use with all patient populations. The product range is suitable for use with adult patients with a body weight as indicated on the individual product pack and the connector of the device. Please refer to the technical data table for size selection. Whilst size selection on a weight basis should be applicable to the majority of patients, individual anatomical variations mean the weight guidance provided should always be considered in conjunction with a clinical assessment of the patient's anatomy.

Warnings

- The i-gel® is supplied in a protective cradle to ensure the device is retained in the correct flexion prior to use and also acts as a base for lubrication. The i-gel® must always be separated from the cradle prior to insertion. The cradle is not an introducer and must never be inserted into the patient's mouth.
- i-gel® must be lubricated according to the instructions for use.
- Remove dentures or removable plates from the mouth before attempting insertion of the device.
- The patient should always be in the 'sniffing-the-morning air' position prior to insertion with the assistant helping to open the patient's mouth, unless head/neck movements are considered inadvisable or are contraindicated.
- Optimum depth of anaesthetic must be achieved prior to attempting insertion (i.e. absence of eyelash reflex, easy up and down movement of the lower jaw, no reaction to pressure applied to both angles of the mandible).
- The leading edge of the i-gel® tip must follow the curvature of the patient's hard palate upon insertion.
- Do not apply excessive force on the device during insertion. It is not necessary to insert fingers or thumbs into the patient's mouth during the process of inserting the device.
- If there is a failure to achieve complete insertion after utilising the standard insertion technique and a jaw thrust, deep rotation or triple manoeuvre has also failed, then the device should be inserted under direct vision by laryngoscopy or one size smaller device should be used.
- After insertion, i-gel® should be taped from maxilla-to-maxilla. The i-gel® should be tied as well as taped during surgical procedures that require the patient to be in a lateral or prone position.
- Excessive air leak during manual ventilation is primarily due to either sub-optimal depth of anaesthesia or sub-optimal depth of i-gel® insertion.
- Particular care should be taken with patients who have an ASA or Mallampati score of III and above or who have fragile and vulnerable dental work, in accordance with recognised airway management practices and techniques.
- As with all supraglottic airways, it is important to ensure the correct size of device is used, lubrication is optimal, the device is inserted and positioned correctly and regularly checked intraoperatively in order to reduce the potential for nerve damage, tongue numbness, cyanosis and other potential complications.
- As with any airway device, when using the i-gel®, it is essential to ensure there is always an adequate expiratory pathway for gas to avoid the potential for barotrauma.
- No attempt should be made to use i-gel® as a conduit for intubation without fiberoptic guidance.
- A gastric tube should not be passed down the gastric channel of the i-gel® if there is an excessive leak through the gastric channel, or in cases of oesophageal trauma, or if there are oesophageal varices or evidence of upper gastro-intestinal bleed, there is a history of upper gastro-intestinal surgery or the patient has bleeding/clotting abnormalities.
- Insertion of a gastric tube in the presence of inadequate levels of anaesthesia can lead to coughing, bucking, excessive salivation, laryngospasm or breath holding.
- Do not attempt to forcibly remove the device if the patient is biting on it. Wait until the patient, on vocal command, has fully opened their mouth or opens their mouth spontaneously.
- Supraglottic airways may be flammable in the presence of lasers and electrocautery equipment.
- The i-gel® incorporates a gastric channel, which when the device is correctly placed, may provide some protection against aspiration of gastric contents. However, the level of protection will vary and cannot be considered definitive.
- In order to reduce the risk of fire, it is essential the user fully complies with the guidance provided by the European Resuscitation Council (ERC), the American Heart Association (AHA) or other appropriate local resuscitation council regarding the safe use of oxygen during defibrillation.
- Whilst size selection on a weight basis should be applicable to the majority of patients, individual anatomical variations mean the weight guidance provided should always be considered in conjunction with a clinical assessment of the patient's anatomy. Patients with cylindrical necks or wide thyroid/cricoid cartilages may require a larger size i-gel® than would normally be recommended on a weight basis. Equally, patients with a broad or stocky neck or smaller thyroid/cricoid cartilage, may require a smaller size i-gel® than would normally be recommended on a weight basis. Patients with central obesity, where the main weight distribution is around the abdomen and hips, might in practice require an i-gel® of a size commensurate with the ideal body weight for their height rather than their actual body weight.
- For patients undergoing Electroconvulsive Therapy (ECT) treatment, a separate bite block specifically designed to be used during ECT must always be used in conjunction with the i-gel®.

Cautions

- Always perform the pre-use checks included in these instructions for use.
- The patency of the device should be checked regularly and particularly after any change in the patients head or neck position.
- i-gel® is a single use device. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection, device malfunction and breakage. Cleaning can lead to device malfunction and breakage.

Size Selection

Select the appropriate size i-gel® by assessing the patient's anatomy and referring to the patient size and patient weight guidance in the Technical data chart below. The i-gel® cuff may look smaller than traditional supraglottic devices with an inflatable cuff of the same numerical size.

Pre-use checks

- Inspect the packaging and ensure it is not damaged prior to opening.
- Inspect the device carefully, check the airway is patent and confirm there are no foreign bodies or a BOLUS of lubricant obstructing the distal opening of the airway or gastric channel.
- Carefully inspect inside the bowl of the device, ensuring surfaces are smooth and intact and also that the gastric channel is patent.
- Discard the device if the airway tube or the body of the device looks abnormal or deformed.
- Check the 15 mm connector fits the patient connection.

Preparation

1. Always wear gloves.
2. Open the i-gel® package and take out the protective cradle containing the device (**figure 1**).
3. In the final minute of pre-oxygenation, remove the i-gel® and transfer it to the palm of the same hand that is holding the protective cradle, supporting the device between the thumb and index finger (**figure 2**). Place a small bolus of a water-based lubricant, such as K Y Jelly, onto the middle of the smooth surface of the cradle in preparation for lubrication (**figure 3**). Do not use silicone based lubricants.
4. Grasp the i-gel® with the opposite (free) hand along the integral bite block and lubricate the back, sides and front of the cuff with a thin layer of lubricant. This process may be repeated if lubrication is not adequate, but after lubrication has been completed, check that no BOLUS of lubricant remains in the bowl of the cuff or elsewhere on the device. Avoid touching the cuff of the device with your hands (**figures 4, 5, 6 and 7**).
5. Place the i-gel® back into the cradle in preparation for insertion (**figure 8**).

Recommended insertion technique

A proficient user can achieve insertion of the i-gel® in less than 5 seconds.

1. Remove the i-gel® from the protective cradle and grasp the lubricated device firmly along the integral bite block. Position the device so that the i-gel® cuff outlet is facing towards the chin of the patient (**figure 9**).
2. The patient should be in the 'sniffing the morning air' position (**figure 9**) with head extended and neck flexed. The chin should be gently pressed down by an assistant before proceeding to insert the i-gel®.
3. Introduce the leading soft tip into the mouth of the patient in a direction towards the hard palate.
4. Glide the device downwards and backwards along the hard palate with a continuous but gentle push until a definitive resistance is felt.
5. At this point the tip of the airway should be located into the upper oesophageal opening (**figure 10a**) and the cuff should be located against the laryngeal framework (**figure 10b**). The incisors should be resting on the integral bite-block (**figure 10c**).
6. i-gel® should be taped down from 'maxilla to maxilla' (**figure 11**).
7. If required, an appropriate size gastric tube may be passed down the gastric channel. See the technical data table below.

Use of the i-gel® as a conduit for intubation

In the known difficult or unanticipated difficult intubation the i-gel® can be used for intubating the patient by passing an appropriate size endotracheal tube (ETT) through the airway channel of the device under fiberoptic guidance. See the technical data table below.

The i-gel® should always be used in accordance with recognised airway management practice for supraglottic airway devices.

Removal

- Once consciousness is regained and protective reflexes such as coughing and swallowing have returned, gently suction around the airway device in the pharynx and hypopharynx. Once the patient is awake or easily arousable with vocal commands, the i-gel® can safely be removed by asking the patient to open his/ her mouth.
- In patients with the possibility of a heightened gag reflex (i.e. smokers, asthmatics or patients with COPD), i-gel® should be removed in deeper planes of anaesthesia and, after removal, the airway maintained with a Guedel airway until protective reflexes have returned and the patient becomes arousable.

Storage conditions

Store in a dry place out of direct sunlight. Products should be stored between -18 °C - +40 °C.

Disposal

Following use, the product must be disposed of in accordance with local hospital infection control and waste disposal regulations.

Technical data

Product code	Size	Patient weight (kg)	Internal volume (ml)	Nominal length of internal pathway (mm)	Maximum size of ET tube that can be inserted down the airway channel (mm)	Maximum size of gastric tube that can be passed down the gastric channel (FG)	Minimum interdentary gap for insertion (mm)	Pressure drop
8203000	3	30–60	22.6	189	6.0	12	27.7	< 0.3cmH ₂ O @ 60L/min
8204000	4	50–90	26.5	194	7.0	12	30.7	< 0.2cmH ₂ O @ 60L/min
8205000	5	90+	34.0	206	8.0	14	32.9	< 0.2cmH ₂ O @ 60L/min

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* Explanation of packing symbols can be found on the Intersurgical website <http://www.intersurgical.com/support> or in standard ISO 15223-1:2016.

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INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Ireland T: +353 (0)1 697 8540 info@intersurgical.ie	España T: +34 91 665 73 15 info@intersurgical-es.com	Nederland T: +31 (0)413 243 860 info@intersurgical.nl	Sverige T: +46 08 514 30 600 info@intersurgical.se	Россия T: +7 495 771 6809 info@intersurgical.ru	South Africa T: +27 (0)11 444 7968 info@intersurgical.co.za	台灣 T: +886 4 2380 5430 info@intersurgical.com.tw	Canada T: +1 905-319-6500 info@intersurgical.ca
France T: +33 (0)1 48 76 72 30 info@intersurgical.fr	Portugal T: +351 219 108 550 info@intersurgical.pt	België/Belgique T: +32 5 2456991 info@intersurgical.be	Danmark T: +45 28 69 80 00 info@intersurgical.dk	Česká Republika T: +420 272 911 814 info@intersurgical.cz	中國 T: +86 519 69817000 info@intersurgical-cn.com	Philippines T: +63 2820 4124 info@intersurgical.ph	Columbia T: +57 1 742 5161 info@intersurgical.com
Deutschland T: +49 (0)2241 25690 info@intersurgical.de	Italia T: +39 0535 20836 info@intersurgical.it	Luxembourg T: +32 5 2456991 info@intersurgical.be	Lietuva T: +370 387 66611 info@intersurgical.lt	Türkiye T: +90 216 468 88 28 info@intersurgical.com.tr	日本 T: 81 (0)3 6863 4359 info@intersurgical.co.jp	USA T: +1 800 828 9633 support@intersurgicalinc.com	Australia T: +61 (02) 8048 3300 info@intersurgical.com.au

patent: <http://us.intersurgical.com/support>

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

9278 Issue 15 RO 07.21