







Laryngeal Tube LTS-D

2nd generation supraglottic airway with drain tube

USE

English

INTENDED USE

The LTS-D is a supraglottic airway with drain tube.

CLINICAL BENEFIT

The LTS-D allows the ventilation of patients.

PATIENT TARGET GROUP

Newborn to adult.

INDICATIONS

- As alternative adjunct to temporarily secure the airway during difficult airway management and cardiopulmonary resuscitation in clinical as well as preclinical environments (current guidelines of the European Resuscitation Council, American Heart Association and American Society of Anesthesiologists).
- For mechanical or spontaneous ventilation of patients during short anaesthesia procedures.

Further indications are not known.

CONTRAINDICATIONS

- Complete obstruction of the upper airway.
- In elective use: patients who have not fasted as well as patients with retained gastric contents.

Further contraindications are not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel with adequate skills in the use of the product.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before use, the products must be inspected visually for damages and a functional check must be performed (see "Visual check / Functional check").
- During MRI treatment, the pilot balloon of the tube can cause artefacts around the area where it is positioned. A non-clinical test verified that the tube can be used in an MRI environment.
- The product is suitable for conducting oxygen, air and anaesthesia gases. The anaesthesia gases used must not interact with the product materials (see material data).
- Laser or electrocaustic equipment may cause the product to be flammable.
- Use original VBM parts only. Using parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.
- The tube and the syringe are for single use and must not be reprocessed. Reprocessing leads to impairment of the function of the products. Any reuse entails a potential risk of cross-contamination.

STERILE EO

- The tube is sterile (ethylene oxide).
- The product must not be used if the packaging is damaged or if shelf life has elapsed.

PRODUCT DESCRIPTION

(a) ISO connector 15 mm O.D.	colour-coded, for easy identification of different sizes
(b) Inflation line	inflates and deflates both cuffs
(c) Pilot balloon	with check valve and luer connection
(d) Valve opener	-
(e) Ventilation	-
(f) Pressure release	-
(g) Distal cuff	seals oesophagus and reduces possibility of gastric insufflation
(h) Ventilation orifices	lie in front of the larynx
(i) Proximal cuff	stabilises the tube and seals naso- and oropharynx
(j) Teeth marks	for easy orientation and repositioning
(k) Drain tube	for gastric tube

USE

- ▶ Choose the appropriate tube size (see table "Product specifications").
- ▶ Remove the valve opener.

Perform visual check and functional check:

- ▶ Check the tube for damage and loose particles.
- ▶ Check the inside of the tube for sufficient patency.

Check for leaks:

- ▶ Deflate the cuffs completely with the syringe (figure 1). A vacuum is created in the cuffs and must be maintained.
- ▶ Alternative: Inflate the cuffs with the syringe.

No air may leak from the cuffs.

Damaged products must not be used.

Tube insertion:

Have a spare tube ready and prepared for immediate use.

- ▶ Deflate the cuffs completely with the syringe (figure 1).
- ▶ Lubricate the tip and posterior aspect of the tube with a water-based lubricant.
- ▶ Pre-oxygenation is advisable in all cases.
- ▶ Head position.

The ideal head position for insertion is the "sniffing position". However, the angle and the shortness of the tube also allows it to be inserted with the head in a neutral position.

CAUTION

Achieve an appropriate depth of anaesthesia before tube insertion!

- ▶ Hold the tube like a pen in the area with the black lines (teeth marks).
- ▶ Open the mouth with the non-dominant hand and apply chin lift (figure 2). The base of the tongue is lifted from the posterior pharyngeal wall allowing a maximum pharyngeal space. Alternatively a jaw thrust maneuver can be used (figure 3).
- ▶ The tongue must not fold back during the downward movement of the tube.

CAUTION

Do not damage the cuffs during tube insertion. The cuffs must therefore be completely deflated.

- ▶ Never use force to insert the tube!

- ▶ Insert the tube either in the midline of the mouth or into the corner of the mouth while rotating the tube.

- ▶ As the proximal cuff passes the incisors, rotate back the tube towards the nose (figure 4) and advance the tube along the base of the tongue until the upper teeth mark (figure 5).
- ▶ Use a syringe to inflate the cuffs with the recommended volume (see table "Product specifications") (figure 6).

Thanks to the high-volume low-pressure cuffs, a seal can already be achieved with low cuff pressures (< 60 cmH₂O).

Ventilation:

- ▶ Connect the tube to a ventilation system (figure 7).
- ▶ Use an appropriate procedure to check ventilation and exclude gastric insufflation.

Gastric tube insertion:

- ▶ To control correct positioning it is recommended to insert a gastric tube (see table "Product specifications") via the drain tube (figure 8).
- ▶ For pressure relief gastric fluids and gases can be aspirated with a gastric tube.

If the gastric tube cannot be advanced this may indicate that the tip of the LTS-D is kinked or a malposition in the piriform fossa.

- ▶ In this case, remove the gastric tube, completely deflate the cuffs and reinsert the LTS-D.

Tips for optimization in case of inadequate ventilation:

1. Tube inserted too deep: Withdraw the tube with inflated cuffs between the teeth marks until ventilation is easy and free flowing (large tidal volume with minimal airway pressure).
2. Tube inserted too shallow: In this case remove the gastric tube, deflate the cuffs completely and reinsert the LTS-D.

Tube fixation:

- ▶ Secure the tube to prevent dislocation.

Cuff pressure measuring:

CAUTION

Excessive cuff pressure can lead to tongue swelling. Colour differences between tongue and lip are an alarm sign for a possible perfusion problem of the tongue:

- Every patient shifting may change the cuff pressure.
- The cuff pressure may change when a cervical spine collar is applied.
- Oxygen, air or anaesthetic gases can change the cuff volume or cuff pressure.
- ▶ Reduce the cuff pressure to the Airway Leak Pressure and check regularly.

Removal of the LTS-D:

- ▶ Deflate the cuffs completely with the syringe (figure 1).

Once the tube is in the correct position, it is well tolerated until the return of protective reflexes.

- ▶ Remove the tube in an area where suction equipment and the ability for rapid intubation are present.

SHELF LIFE

Expiry date: see product label.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed in accordance with the applicable national and international regulations.

PRODUCT SPECIFICATIONS

Tube size	#0	#1	#2	#2.5	#3	#4	#5
 Patient weight	< 5 kg	5 - 12 kg	12 - 25 kg	-	-	-	-
 Patient height	-	-	-	125 - 150 cm	< 155 cm	155 - 180 cm	> 180 cm
Colour coding ISO connector 15 mm	transparent	white	green	orange	yellow	red	purple
Recommended cuff volume	10 ml	20 ml	35 ml	40 - 45 ml	50 - 60 ml	70 - 80 ml	80 - 90 ml
Maximum cuff pressure	60 cmH ₂ O						
Maximum size gastric tube	10 Fr		16 Fr		18 Fr		
Bronchoscopy via ventilation lumen	< 3.0 mm		< 4.0 mm		< 6.0 mm		
External diameter of the tube	9 mm		14 mm		17.6 mm		

SCOPE OF DELIVERY

REF	Description					Box
	Larygeal Tube LTS-D Box of 1					
32-06-100-1	#0	Newborn	< 5 kg	transparent	colour-coded syringe 20 ml	1 piece
32-06-101-1	#1	Infant	5 - 12 kg	white	colour-coded syringe 20 ml	1 piece
32-06-102-1	#2	Child	12 - 25 kg	green	colour-coded syringe 60 ml	1 piece
32-06-125-1	#2.5	Child	125 - 150 cm	orange	colour-coded syringe 60 ml	1 piece
32-06-103-1	#3	Adult	< 155 cm	yellow	colour-coded syringe 60 ml	1 piece
32-06-104-1	#4	Adult	155 - 180 cm	red	colour-coded syringe 100 ml	1 piece
32-06-105-1	#5	Adult	> 180 cm	purple	colour-coded syringe 100 ml	1 piece
	Larygeal Tube LTS-D Set of 10					
32-06-000-1	#0	Newborn	< 5 kg	transparent		10 pieces
32-06-001-1	#1	Infant	5 - 12 kg	white		10 pieces
32-06-002-1	#2	Child	12 - 25 kg	green		10 pieces
32-06-025-1	#2.5	Child	125 - 150 cm	orange		10 pieces
32-06-003-1	#3	Adult	< 155 cm	yellow		10 pieces
32-06-004-1	#4	Adult	155 - 180 cm	red		10 pieces
32-06-005-1	#5	Adult	> 180 cm	purple		10 pieces
	LTS-D Emergency set					
32-06-309-1	#0 + #1 + #2 + #2.5		Child		colour-coded syringe 60 ml	1 piece each
32-06-209-1	#3 + #4 + #5		Adult		colour-coded syringe 100 ml	1 piece each
Accessories						
54-04-111	Syringe 20 ml				colour-coded	10 pieces
54-04-444	Syringe 60 ml				colour-coded	10 pieces
54-04-555	Syringe 100 ml				colour-coded	10 pieces

MATERIAL DATA

Tube Cuffs Inflation line Pilot balloon	PVC (polyvinyl chloride)
Tube connector	PC (polycarbonate)
Check valve and luer connection	PVC (polyvinyl chloride), PBT (polybutylene terephthalate), NBR (nitrile butadiene rubber), stainless steel
Syringe	PP (polypropylene), IR (isoprene rubber)

Conversion of pressure units: 1 hPa = 1.01973 cmH₂O = 0.75006 mmHg
Conversion of length units: 1 Fr = 0.33 mm

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