Instructions for use

PerSys Medical

Instructions for use


PerSys Medical

Instructions for use

PerSys Medical

Instructions for use

PerSys Medical

Instructions for use

PerSys Medical

Instructions for use


PerSys Medical

Instructions for use

PerSys Medical

Instructions for use

PerSys Medical

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ENGLISH

NIO ADULT DESCRIPTION

The NIO Adult is a single-use, automatic, spring-loaded intraosseous sterile device with dual safety features for enhanced patient and caregiver safety. After deployment, its built-in stabilizer base, in conjunction with a stand-alone NIO Fixation dressing, ensures that the needle is firmly secured.

INDICATIONS FOR USE/INTENDED USE

The NIO Adult is intended to provide intraosseous (IO) access in the proximal tibia and the humeral head as an alternative to emergency IV access. Per urgent medical necessity, humeral head IO access may be used when rapid fluid or pharmacological resuscitation is required. The device is intended for use in adult patients only (above 12 years of age).

NIO ADULT INSTRUCTIONS FOR USE

Always wear gloves during the intraosseous access procedure. Disinfect the skin at the insertion site, per institutional protocol.

1. Open the blister pack by peeling off the Tyvek and take out the NIO Adult (NIO-A). Make sure that the NIO-A is free of all packaging parts.

2. Select one of the following insertion sites:

2A. Primary site for intraosseous insertion: proximal tibia

Approximately 2cm (1in) medial and 1cm (0.5in) proximal to the tibial tuberosity.

2B. Secondary site for intraosseous insertion: humeral head

Adduct the patient's arm and locate the greater tubercle, next to the head of the humerus. Palpate the surgical neck of the humerus with your thumbs. Move 1cm (0.5in) towards the greater tubercle. Note: To prevent accidental removal of the device, immobilize the patient's limb throughout the procedure.

3. **Note:** This is a two-handed procedure. With your non-dominant hand, hold the NIO-A by the textured dots and position the device at the insertion site, at a 90° angle to the skin. Your non-dominant hand should maintain this position and continue to **stabilize** the device for the duration of the procedure.

4. **Unlock** the NIO-A by rotating the cap 90° in either direction.

5. Place the palm of your dominant hand over the cap and apply downward pressure, **pressing** the device against the patient's skin. While maintaining downward pressure, **squeeze** the trigger wings upward. This will deploy the device.

6. With your non-dominant hand, hold the stabilizer base against the patient's skin. With your dominant hand, gently lift the NIO-A in an upward and rotating motion. At this point, it is recommended to use the NIO Fixation dressing to affix the stabilizer base to the patient's skin. Remove the paper liner from the underside of the NIO Fixation dressing to expose its adhesive tabs. Use your fingers to affix the adhesive tabs to the patient's skin, ensuring that the NIO Fixation dressing is firmly adhered.

7. While holding the blue stabilizer base and cannula in place, remove the trocar by pulling it up (a twisting motion may be necessary). The trocar removal notch, located at the distal end of the NIO-A, can assist in removing the trocar from the cannula. Once removed, place the trocar into an appropriate biohazard container.

8. Connect a syringe and confirm a secure fitting. If desired, aspirate for bone marrow. Always confirm successful needle placement by flushing with up to 20cc of fluid, or per institutional protocol. Failure to properly flush the cannula may result in limited or no flow. Connect any standard system for infusion.

Note: Be sure to steady the stabilizer base when connecting/disconnecting any Luer lock line or syringe to/from the cannula. Do not leave the hub open to air.

9. Removal instructions: While holding the stabilizer base, loosen the NIO Fixation dressing by pulling its adhesive tabs away from the skin. Remove the stabilizer base by twisting and pulling in an upward motion. Dispose of the removed components in an appropriate biohazard container. Hold manual pressure until hemostasis is achieved. Cover the insertion site with a sterile occlusive dressing. On the occlusive dressing, write the time and date of removal, as well as "NO IO x48hrs" (or follow your institutional protocol).

CONTRAINDICATIONS

Stop use of the device should any of the following diagnoses or occurrences arise: • Fracture of the bone selected for insertion or fracture of any bone within the limb selected for insertion. • Infection at the insertion site. • Previous orthopedic procedure near the insertion site, prosthetic limb or joint. • Any previous intraosseous attempt in the same bone within 48 hours. • Inability to locate anatomical landmarks or the presence of excessive tissue at the insertion site. • Tumor at the insertion site. • Osgood-Schlatter disease. • Deformation of the insertion site. • Abnormalities of bone strength (e.g., osteogenesis imperfecta, osteopetrosis, osteoporosis).

WARNINGS

• Read this Instructions for Use leaflet carefully • Do not use if package is damaged. • Use aseptic technique. • Discontinue infusion if any signs of infiltration are apparent, including tissue swelling around the insertion site. • Do not position the NIO-A toward the joint space or epiphyseal plate. • Use of the NIO-A is restricted to skilled and authorized medics, nurses, paramedics, and doctors who are trained to use the device. • Do not reuse or re-sterilize the device. Reusing the disposable device may lead to infection, mechanical failure, or harm to the operator and/or patient. • This device has a limited shelf life. Do not use beyond the expiration date printed on the device label. Using an expired device may lead to device malfunction or patient injury. • When using any intraosseous device, the possibility of air embolism, fat embolism, fracture of the targeted bone, infection, and extravasation exists. • The NIO-A contains exposed sharp components that should be disposed of in an appropriate medical biohazard waste container. Note: Dispose of the used device in accordance with applicable local laws and regulations, per institutional protocol. • The NIO-A consists of a metal needle and is not MRI compatible. • The NIO-A consists of a metal needle, which may cause scatter artifacts on computed tomography (CT) scans. • The safe use of the NIO-A in patients with osteoporosis, osteopetrosis, Osgood-Schlatter disease, or other bone pathology or deformities has not been proven. These conditions may obscure landmarks. • Caution: This device is to be used on the order of, and by, a licensed physician or licensed practitioner. Rx only.

RECOMMENDATIONS

• Recommended penetration depth for proximal tibia or humeral head 2.5cm (1in). • Reassess the IO site immediately following step 8 (confirm needle placement). After beginning drug administration, it is recommended to frequently monitor the limb, every ten minutes for at least the first half-hour. • Use a pressure bag for optimal infusion rates. • Pre- and post-infusion, regularly monitor the limb for complications. Prior to drug administration, the cannula should be verified for placement and patency by confirming its stability in the bone. Bedside Doppler ultrasonography may be utilized to ascertain IO placement and flow. • For conscious patients, consider a local anesthetic, such as lidocaine, per institutional protocol. • To prevent accidental removal of the IO following insertion into the humeral head, immobilize the patient's arm. • The cannula should be removed within 24 hours of insertion or once alternative venous access has been established.

RECOMMENDED STORAGE CONDITIONS

Store at room temperature. The device should be handled with care.

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