LOW PRESSURE REGULATORS; INSTRUCTIONS FOR USE
NÍZKOTLAKÉ REDUKČNÍ VENTILY; NÁVOD K POUŽITÍ
LES DÉTENDEURS À BASSE PRESSION; NOTICE D'UTILISATION
LAGEDRUKREGELAARS; HANDLEIDING
AZ ALACSONYNYOMÁSÚ NYOMÁSCSÖKKENTŐ SZELEPEK; KEZELÉSI ÚTMUTATÓ
NIEDERDRUCK-DRUCKMINDERVENTILE; GEBRAUCHSANLEITUNG
REDUKTORY NISKOCIŚNIENIOWE; INSTRUKCJA UŻYCIA
НИЗКОНАПОРНЫЕ РЕДУКЦИОННЫЕ КЛАПАНЫ; ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LOS REGULADORES DE BAJA PRESIÓN; INSTRUCCIONES DE USO
OS REGULADORES DE BAIXA PRESSÃO; MANUAL DE UTILIZAÇÃO
I REGOLATORI DI BASSA PRESSIONE; ISTRUZIONI PER L'USO
LÅGTRYCKSRREGULATORERNA; BRUKSANVISNING
LAVTRYKKSREDUKSJONSVENTILER; BRUKSANVISNING
LAVTRYKSREDUKTIONSVENTILER; BRUGSANVISNING
MATALAPAINISET PAINEENSÄÄTIMET; KÄYTTÖOHJE
ZEMSPIEDIENA REGULATORI; LIETOŠANAS PAMĀCĪBA
ŽEMO SLĖGIO REDUKCINIAI VOŽTUVAI; NAUDOJIMO INSTRUKCIJA
MADALRÖHUREGULAATORID; KASUTUSJUHEND
REDUCTOARE DE JOASĂ PRESIUNE; INSTRUCȚIUNI DE UTILIZARE
NÍZKOTLAKOVÉ REDUKČNÉ VENTILY; NÁVOD NA POUŽITIE
NISKOTLAČNI REDUKCIJSKI VENTILI; UPUTE ZA UPORABU
NIZKOTLAČNI REDUKCIJSKI VENTILI; NAVODILA ZA UPORABO
ΡΥΘΜΙΣΤΕΣ ΧΑΜΗΛΗΣ ΠΙΕΣΗΣ; ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
РЕДУЦИР ВЕНТИЛИТЕ ЗА НИСКО НАЛЯГАНЕ; ИНСТРУКЦИЯ ЗА ЕКСПЛОАТАЦИЯ
低压调压器; 使用说明
DÜŞÜK BASINÇ REGÜLATÖRLERİ; KULLANIM TALİMATLARI
ادخست سال تامی ل عمل ال؛ اضیفی نزل طغض ال؛ اضیفی تامام ص

1. FOREWORD

GCE Medical Regulators are medical devices classified as class IIb according to the Medical Device Regulation (EU) 2017/745.

Compliance with the General Safety and Performance Requirements of (EU) 2017/745 Medical Device Regulation is partly based upon EN ISO 10524-4 standard.

Note: *The precautions/warnings specified for oxygen are valid for all gases (gas specification is stated in Appendix Nr. 1), if not specified differently.*

2. INTENDED USE

Low pressure regulators (Mediflow® Ultra II) are designed to be connected directly at terminal units of a medical gas pipeline system or gas-specific connection points of devices such as medical high pressure regulators. Low pressure regulators fitted with the flow-selector device are used to select specific outlet flow suitable for administering medical gases in the treatment of patients.

Low pressure regulators are designed to deal with fluctuating inlet pressure which causes the medical device to deliver accurate and stable flow consistently.

Low pressure regulators are not intended to be used with gases for driving surgical tools.

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

-  Always maintain oxygen cleanliness standards:
 -  Keep the medical device and its associated equipment away from:
 - heat sources (fire, cigarettes, ...),
 - flammable materials,
 - oil or grease (take a great care in the use of hand creams),
 - water,
 - dust.
 -  The medical device and its associated equipment must be prevented from falling over.
 -  Use only the medical device and its associated equipment in well ventilated area.
- Before initial use the medical device should be kept in its original packaging.

GCE recommends use of the original packaging (including internal sealing bag and caps) if the medical device is withdrawn from operation (for transport, storage).

Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20/+60 °C		-30/+60 °C

 In case of storage at a temperature below -20 °C, do not operate the medical device until it has been allowed to increase its temperature to a minimum of -20 °C.

Note: Low pressure regulators do not indicate an occurring flow or its accuracy. The accuracy of the outlet flow may be influenced by varying outlet resistance, or variations in the ambient temperature from 0°C to +40°C.

4. PERSONNEL INSTRUCTIONS

The medical device provider must ensure that all personnel handling the medical device are provided with the operating instructions & performance data. Gas source is not an integrated part of the medical device.

4.1. INTENDED USER

The intended user is any person, i. e. Healthcare professional or lay user*, who has read and understood the instructions for use and is able to operate the device correctly.

*The device is intended for lay users whose ability to operate the device is approved by the healthcare professional who prescribed the gas therapy.

 Do not use the medical device without proper familiarization with the medical device and its safe operation as defined in this Instructions for use. Be aware of particular information and knowledge required for the gas in use.

Note: In case of doubt regarding the use of the device as described by this Instructions for use, contact the medical device provider or manufacturer.

Note: In case of doubt regarding the gas therapy, contact your healthcare professional.

5. MEDICAL DEVICE DESCRIPTION

Gas from the medical gas pipeline system enters through the inlet quick connector, is then calibrated in the flow control head and delivered to the user through the flow outlet user connection. Different flows can be selected by the flow control knob on the flow head.

A - INLET QUICK CONNECTOR (QC)

The low pressure regulator is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe.

B - FLOW CONTROL HEAD

Different flow rates can be selected by turning according to the “+/-” marking on the flow control knob, which is part of the flow head. The flow control head is fitted with a stop function in between maximum flow and zero position. It is not possible to directly set the flow from the zero position to the maximum flow.

C - USER OUTLET

The user outlet (flow outlet) connection is the outlet from the flow control head. It is fitted with a fir-tree (for hose nipple) or threaded (for humidifier) connector. This outlet is intended for supplying a gas flow (l/min) at atmospheric pressure directly to a patient, e.g. through a cannula or a facemask.

 **The flow outlet shall not be used for driving any medical device.**

Note also that medical devices colour (specially flow control knob) might not follow any gas colour coding.

Note: Detailed information about the performance and technical data related to your medical device can be found in Appendix Nr. 1. Appendix Nr. 1 is in form of symbols. Explanations for the symbols are available in Chapter 12 - Glossary.

EN

6. INSTALLATION

Not applicable for this medical device.

7. OPERATIONS

 The medical gas is not an integrated part of the medical device. Gas therapy must be prescribed by a healthcare professional only.

7.1. BEFORE USE

EN

VISUAL INSPECTION BEFORE USE

- Check if there is visible external damage to the medical device (including medical device labels and marking) and on the gas cylinder. If it shows signs of external damages, remove it from service and identify its status.
- Remove caps from inlet and/or flow outlet. Keep caps in a safe place for reuse during transport or storage.
- Visually check if the medical device or the medical gas cylinder is contaminated; if needed, for the regulator, use the cleaning procedure detailed in this section (if required for the cylinder, refer to the gas cylinder manufacturer cleaning recommendation).
- Check that the total lifetime of the medical device and the gas cylinder has not been exceeded, (refer to GCE or owner's date coding system). If lifetime has been exceeded, remove the medical device (or the gas cylinder) from service and suitably identify its status.
- Ensure that the medical device inlet stem is compatible with the medical cylinder valve (gas/ thread type).
- Set flow control knob on ZERO position - until the flow control knob engages.

 The medical device is dedicated only for use with the gas specified on its labelling. Never try to use for another gas.

FITTING TO MEDICAL GAS PIPELINE SYSTEM TERMINAL UNITS

Low Pressure Regulators shall be connected in accordance with a suitable type of terminal unit or with a regulator with an appropriate quick connector, see appendix Nr. 2.

LEAKAGE TEST

After connection to medical gas pipeline system terminal unit or to medical regulator, perform visual and audible check for possible leakages:

- medical device inlet connection to terminal unit or to medical regulator,
- inlet connector to medical device body,
- user outlet.

 If any leakage is detected, use the procedure in chapter 7.3 and return the medical device to GCE for service.

FUNCTIONAL TEST

- Check that there is a gas flow at each setting (for instance, by listening for the sound of gas flow or checking presence of bubbles in an humidifier).
- Reset flow control knob on ZERO position.

7.2. CONNECTION & USE OF OUTLET

Associated equipment, see chapter 8.

-  **Before connecting any associated equipment or medical device to the medical device, always check that it is fully compatible with the medical device connection features & the medical device performances.**

FLOW OUTLET CONNECTION

-  **When connecting any associated equipment to the flow outlet make sure that it is not connected to the patient before operating the medical device.**
- Ensure the hose/humidifier is compatible with the flow outlet feature.
 - Push the hose onto the regulator flow outlet/outlet for humidifier.
 - Ensure the hose/humidifier is well engaged.

FLOW SETTING

- Set the medical gas flow control knob on one of the available flow rates by turning the flow control knob in the requested position.
-  **Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position or in zero position.**
-  **Always ensure that the flow control knob has engaged and not placed between two settings otherwise the medical device will not deliver the correct flow of medical gas.**
-  **During use ensure that the connected hose is not mechanically damaged, throttled, or broken, which would reduce or break the supply of gas to the patient.**

7.3. AFTER USE

- Before disconnecting the low pressure regulator from the medical pipe line system terminal unit or medical regulator, reset the flow control knob to ZERO position.
- Vent gas pressure from downstream equipment (i.e. via the connected consumers).
- Remove connections from user outlets.
- Refit pressure outlet and flow outlet protection caps. Before refitting the caps, ensure they are clean.

8. ASSOCIATED EQUIPMENT

To be connected to the flow outlet:

- Humidifiers, breathing masks or cannulas, gas savers, nebulizers.

To be connected to the inlet:

- Terminal units of medical gas pipeline systems complying with ISO 7396-1 or terminal unit attached to pressure regulator complying with ISO 10524-1 or ISO 10524-3.

9. CLEANING

The medical device is not sterile and does not require cleaning for contamination control. The purpose of cleaning is solely to remove contamination on the external surfaces (e.g. physical dirt) that could potentially cause safety issues with the medical device.

EN

Remove dirt with a soft cloth dampened in oil free soap water & rinsed with clean water. Cleaning can be carried out with an alcohol-based solution (with dampened wipes). If other cleaning solutions are used, check that they are not abrasive and they are compatible with the medical device materials (including labels) and gas (convenient cleaning solution - i.e. Meliseptol).

Do not use cleaning solutions containing ammonia!

Do not immerse in water or any liquid.

Do not expose to high temperatures (such as an autoclave).

To apply the cleaning solution, do not spray it, as the spray may enter the inner parts of the medical device and cause contamination or damage.

Do not use pressure water as it could damage or contaminate the medical device.

Do not use the medical device under any circumstances if the inner parts of the medical device have been contaminated. It must be withdrawn from service.

10. MAINTENANCE

MAINTENANCE

No special maintenance or service, apart from the tests before use, is needed. However, to make sure the medical device is in good working order, it would be good if the owner/cylinder distributor performs the checks (see 7.1) himself on a regular basis (ex. once every second year) and/or at every cylinder exchange. This is just for the owner to ensure that the medical device works well, especially in cases where the user has some health problems and is not able to check the medical device himself properly.

10.1. REPAIR

Repairs activities cover the replacement of the following damaged or missing components (spare parts):

- inlet and outlet connector,
- piston.

The repairs shall be carried out by a GCE authorised person only.

Any medical device sent back to a GCE-authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair, maintenance). For medical devices to be repaired a short description of the fault and any reference to a claim number might be helpful. Some repair activities concerning the replacement of the damaged or missing components can be carried out by the owner of the medical device.

The following spare parts can be replaced only:

- caps,
- flow knob and stickers,
- hose nipple (including o-ring).

- ⚠ **To use a suitable component, contact our service organization.**
- ⚠ **All labels on the medical device must be kept in good, legible condition by the owner and the user during the entire medical device lifetime.** **EN**
- ⚠ **All seals and O-rings must be kept in a dry, dark, and clean environment by the owner and the user during the entire medical device lifetime.**
- ⚠ **Use only original GCE components!**

10.2. SERVICE

GCE recommends a medical device periodic inspection to be undertaken at least once a year to do leakage test and check the proper functionality of the medical device.

10.2.1. LEAKAGE TEST

- After connection to medical gas pipeline system terminal unit or to a medical regulator, perform visual and audible check for possible leakages:

LEAKAGE TEST OF THE FLOW SELECTOR FOR FLOW SELECTION:

1. Select the value „0“ by turning the handwheel of the flow selector. When closing don't use excessive torque.
2. Connect the flow outlet and hose nipple with a hose, immerse the end of the hose into a jar filled with water.
3. Connect the flow selector to the source of compressed air max. 8 bar.
4. The test: leave the flow selector 1min under pressure. Number of bubbles shall not be higher than 3 bubbles per 1 min.
5. When test ready, disconnect the flow selector from the source of compressed air.

EXTERNAL LEAKAGE TEST:

1. Select the value „0“ by turning the handwheel of the flow selector. When closing don't use excessive torque.
2. Close the outlet connection of the flow selector using a plug, the plug has to have the same thread as the outlet connector of the flow selector or, in cases where the hose nipple is used, use the hose with the blind plug.
3. Connect the flow selector to the source of compressed air max. 8 bar.
4. Select the maximal value by turning the handwheel of the flow selector. When opening don't use excessive torque.
5. Possible leakage can be detected by application of soap water on medical device joints.
6. After the test disconnect the flow selector from the source of compressed air, remove the plug from outlet connection and carefully clean the surface of the flow selector, see chapter 9 (cleaning).

- ⚠ **If any leakage is detected, use the procedure in chapter 7.3 and return the medical device to GCE for service.**

 The medical device shall not be disassembled or repaired when connected to the source of the pressure.

11. MEDICAL DEVICE LIFETIME

11.1. LIFETIME

EN

Maximum lifetime of the medical device is 10 years.

At the end of the medical device's lifetime, the medical device must be withdrawn from service (max. 10 years).

11.2. SERIAL NUMBER AND DATE OF MEDICAL DEVICE

MARKING PRIOR TO IMPLEMENTING UDI CODING

Form of nine digits serial number stamped on the medical device is following: YYMMXXXXX

YY: year of production

MM: month of production

XXXXX: sequence number

Example: serial number 290300521 shows the medical device produced in March 2029, with sequence number 521.

11.3. MARKING UDI CODING IMPLEMENTATION

There are 4 marked identifiers on the medical device surface:

14 numbers with prefix (01): GTIN - Global Trade Identification Number (identification of the selector type)

10 numbers with prefix (10): batch number (LOT)

9 numbers with prefix (21): serial number in format YYMMXXXXX

6 numbers with prefix (17): the end of the medical device's lifetime in format YYMMDD

Example: marking number (17) 290227 shows the medical device with end of life time on 27th of February 2029.

As the integral part of UDI marking is a Data matrix, which consists of the same data.

11.4. DISPOSAL OF WASTE MANAGEMENT

The owner of the medical device shall prevent the reuse of the medical device and handle the medical device in compliance with the "Directive of European Parliament and Council 2008/98/EC on waste".

Contact your local GCE representative for further details, before returning or removing the medical device. All components should be disposed of according to local environmental laws applicable in the country of disposal.

If required, a list of component materials is available from GCE.

11.5. REACH AND RoHS

In accordance with Article 33 of REACH, GCE, s.r.o., as the legal manufacturer, shall inform all customers if any materials exist contain 0.1% or more of substances included in the list of Substances of Very High Concern (SVHC). The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-100-4, CAS no. 7439-92-1. The lead will not be released into the gas or surrounding environment during normal use. After the end of life, the medical device shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact on the environment and health. To date, we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE medical device.

11.6. REPORTING OF A SERIOUS INCIDENT

If any suspected serious incident has occurred in relation to this medical device, notify the manufacturer, e-mail: adverse_events@gcegroup.com, and the competent authority of the State in which the user and/or patient is established.

By reporting a suspected serious incident, you can help obtain more information about the safety of this medical device.

12. GLOSSARY

	Consult instructions for use		Suitable for Hospital care use
	Caution		Suitable for Home care use
	Keep away from heat and flammable material		Suitable for Emergency care use
	Keep away from oil and grease		Serial number
	Keep dry		Catalogue number
	Use by date		Customer number
	Inlet pressure range		Batch code
	Outlet pressure		Fragile, handle with care

	Max outlet pressure (closing pressure)		Manufacturer
	Inlet parameter		Weight of medical device
	Temperature limit		Outlet parameter
	Contains hazardous substances		Authorised representative for Switzerland
	Medical Device		Unique Device Identification

13. WARRANTY

The Standard Warranty period is 2 years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the medical device manufacture shown on the medical device).

The standard warranty is only valid for medical devices handled according to Instructions for use (IFU) and general industry good practice and standards.

APPENDIX:

Nr. 1 - Technical and performance data

Nr. 2 - Quick coupling feature and connecting / disconnecting procedure

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© GCE, s.r.o.