



# ri-pen®

## Diagnostic pupil penlight

Sustainable penlight for medical diagnosis with LED illumination

- Photobiologically validated for lighting of pupils according to EN 62471:2008
- UV free light through LED illumination
- Better vision due to a centered light beam
- Nickel free
- Improved hygiene – desinfectable
- Light aluminium casing with metal clip
- Simple on/off function
- Including batteries
- 2-year warranty



### Order information:

#### ri-pen®, pack of 1

Art. Nr.: 5075-526 black  
Art. Nr.: 5074-526 silver  
Art. Nr.: 5077-526 red  
Art. Nr.: 5071-526 blue  
Art. Nr.: 5078-526 purple  
Art. Nr.: 5072-526 green

#### ri-pen®, pack of 6

Art. Nr.: 5070 assorted colours  
Art. Nr.: 5075 black  
Art. Nr.: 5074 silver  
Art. Nr.: 5077 red  
Art. Nr.: 5071 blue  
Art. Nr.: 5078 purple  
Art. Nr.: 5072 green

### Technical data:

3 V LED  
5.300 Kelvin  
40.000 Lux max.  
Field of view approx. 40 mm diameter  
(distance 15 cm)

**Batteries:** 2 x Alkaline AAA, 1,5 V